## MTC Bio Culture & Centrifuge Tubes



- Round Bottom Culture Tubes
- Conical Centrifuge Tubes
- High Clarity Polymers
- Printed Graduation & Marking Spot
- Sterile; RNase & DNase Free



These disposable sterile tubes are made of high quality polypropylene (PP) that will withstand temperatures from -190° to +120°C. They are also available in crystal clear polystyrene (PS) with a more limited range of 0° to 70°C. Caps are made of polyethylene (PE) and have the popular two position design (up/ loose for aerobic work and down/ sealed for anaerobic cultures). Polystyrene tubes may be centrifuged at up to 2000xg. Polypropylene tubes withstand higher forces, up to 4000xg.

## Sterile Centrifuge Tubes

Conical 15ml and 50ml centrifuge tubes are made of premium polypropylene for clarity and centrifugal strength. Tubes and caps may be frozen at temperatures as low as -86°C. The maximum centrifugal force is 12,000xg for the 15ml tubes iand 15,000xg for the 50 ml tubes\*. Caps have a flat top surface for easy marking and labeling. Cat. No. C2605 features a "skirted" bottom that allows the tube to stand upright without a tube rack.

\*Refers to liquids w/ a specific gravity of 1.0, using appropriate rotors, cushions & adapters.



Ph: 908-279-0280 www.mtcbio.com

## **Ordering Information**

<b>Culture Tubes</b> (supplied with dual position caps attached)										
Cat. No	Capacity	Dimensions	Mate- <u>rial</u>	Packaging Unit	<u>Units</u> Per Case	Price Per Case				
T8301	5ml*	12x75mm	PP	25 in heat seal bag	40	\$109.00				
T8305	5ml*	12x75mm	PS	25 in heat seal bag	40	\$109.00				
T8320	5ml*	12x75mm	PP	125 in heat seal bag	8	\$99.00				
T8300	5ml*	12x75mm	PS	125 in heat seal bag	8	\$99.00				
T8235	14ml*	17x100mm	PP	25 in notched bag	40	\$209.00				
T8260	14ml*	17x100mm	PP	1 per sterile peel	500	\$129.00				
T8225	14ml*	17x100mm	PS	25 in notched bag	40	\$215.00				
T8250	14ml*	17x100mm	PP	1 per sterile peel	500	\$129.00				

<u>Centrifuge Tubes</u> (supplied with screw caps attached)

Cat. No	Capacity	Dimensions	Mate- rial	Packaging Unit	<u>Units</u> Per Case	Price Per Case			
C2600	15ml	17x118mm	PP	50 tubes in foam rack	10	\$129.00			
C2601	15ml	17x118mm	PP	50 tubes in notched bag	g 10	\$119.00			
C2602	50ml	29x115mm	PP	25 tubes in foam rack	20	\$149.00			
C2603	50ml	29x115mm	PP	25 tubes in notched bag	g 20	\$139.00			
C2603-F	3 50ml	29x115mm	PP	25 tubes in notched bag	g 20	\$139.00			
C2605*	50ml	29x115mm	PP	25 tubes in notched ba	g 20	\$159.00			
Tube Racks (for all standard size tubes)									
R1070 Rack for 30x15ml tubes, pack of 5									
R1030	$eq:control_co$								

<sup>\*</sup> No. 2605 is a self-standing, sterile tube