## Sureseal S





Syringe Permeable
2mm circle is thinner,
allowing through-cap injection

Raised locking rim seals into groove in sidewall

Smooth Rounded Lip prevents "red thumb syndrome"

Glass-clear side walls maximum transparency BioClear™ polypropylene

Ultra-smooth interior surface Nano-polished molding process

Advanced *BioClear*™ polypropylene enhances pellet vizualization and minimizes optical distortion

Ring and Groove closure protruding ring on cap locks into groove in upper sidewall with a distinct "click"

Frosted marking surface on non-graduated side

Vivid graduations at 0.1ml and every 0.25ml

Reinforced Conical Bottom withstands extreme g-forces

- Sterile comes <u>Standard</u>
   by electron beam radiation
- Self-Standing bag with gusseted bottom
- Stop-Pops<sup>™</sup> Supplied
   Stop-Pop cap locks in every bag







Also available in these great colors!

MTC Bio's **SureSeal S<sup>TM</sup>** microcentrifuge tubes are designed for a broad range of research applications, including enzyme digests, DNA isolation/purification, centrifugation, incubation and general sample/reagent storage. With a variety of design advances, **SureSeal** Tubes improve sample visualization, safety, stability and overall ease of use.

SureSeal S tubes are manufactured under clean room conditions, using only BioClear™ co-polymer, an advanced polypropylene that maximizes clarity. BioClear has been formulated to be crystal clear, without compromising centrifugal strength (tested up to 28,000xg) or the usable temperature range (-90°C to +121°C).

The oversized tube cap (with frosted marking surface), combined with the unique ring & groove closure, enhances the security of the liquid-tight seal, while minimizing the required locking forces for a more ergonomic opening/closing of the tube caps.

2600A Main Street Extension Sayreville, NJ 08872 Ph: 908-279-0280 www.mtcbio.com

All SureSeal S tubes are RNase, DNase and pyrogen free.

## Ordering Information

Cat. No.	<u>Description</u>	<u>Qty</u>
C2000*	1.5mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
C2000-AST	1.5mL microtube, sterile, self-standing bag, 10 Stop-Pops, assorted colors	500
C2002	2.0mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
C2007	0.5mL microtube, clear, sterile, self-standing bag,	500
C2007-AST	0.5mL microtube, sterile, self-standing bag, assorted colors	500

\* Please add (-B) for blue (-R) for red,(-G) for green, (-O) for orange, (-OB) for opaque black, (-P) for purple and (-Y) for yellow