

# ProPette™ MPC

## MANUAL PIPETTE CONTROLLER



- **For glass or plastic pipettes**
- **Accommodates pipettes up to 100ml**
- **Easy, single lever control**
- **Autoclavable lower body**
- **Fingertip speed control**

ProPette™ MPC is the precise and easy to use manual pipette controller that's always ready to use. After attaching a pipette, one squeeze of its large rubber filling bulb is all it takes to start pipetting - there are no wires, no tubing and no battery to recharge! Begin aspirating by touching the EZ-Flow™ lever which is designed to precisely control the filling speed, based on how gently or firmly you press the lever.

A built-in replaceable filter ensures optimum flow rate and protects the controller from liquid penetration and over-fill. A separate button provides a positive blow-out function when required.

Made to accept pipettes from 0.1 to 100ml, the silicone "chuck" adapter is flexible, completely autoclavable, and will accommodate all types of serological, "blow-out", "to deliver," glass and plastic pipettes with a soft, precise, air-tight seal.

**mtc-Bio**

2600A Main Street Extension  
Sayreville, NJ 08872

Ph: 908-279-0280 [www.mtcbio.com](http://www.mtcbio.com)

### Ordering Information

| <u>Cat. No.</u> | <u>Description</u>                                  | <u>Price</u> |
|-----------------|---|--------------|
| P6420           | ProPette MPC, Manual Pipette Controller with filter | \$69.00      |
| P6080-FI        | Replacement filter, 0.45um, package of 5            | \$20.00      |