

# Individually Wrapped & Sterile 15 & 50mL Conical Tubes



- **Sterile, and individually wrapped**
- **RNase, DNase, and Pyrogen free**
- **Easy peel wrappers**
- **Flat markable screw-caps**

MTC Bio has taken the same screw-cap conical tubes that you know and trust and repackaged them in to individual sterile packets. No more fumbling around with open bags or racks of 50 tubes. No more wondering if the tubes in the bag you opened yesterday are still sterile. With the individual wrapping, each tube is sterile right up until the moment you open the convenient peel pack.

Like our standard bagged and racked tubes, the individually wrapped conical tubes are made of premium polypropylene for clarity and centrifugal strength. Tubes and caps may be frozen at temperatures as low as  $-86^{\circ}\text{C}$ . The maximum centrifugal force for the 15mL tubes is 12,000xg, while the 50mL tubes can reach 15,000xg. Caps have a flat top surface for easy marking and labeling. Individual tubes are sterile, DNase, RNase, and pyrogen free, and include high contrast graduations, and a white "frosted" marking area.

**MTC Bio**

2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 [www.mtcbio.com](http://www.mtcbio.com)

## Ordering Information

<u>Cat. No</u>	<u>Description</u>	<u>Price</u>
C2615-W	15mL centrifuge tubes, individually wrapped, sterile, 300/pk	\$79.50
C2650-W	50mL centrifuge tubes, individually wrapped, sterile, 300/pk	\$99.50