## **SUMMER / FALL** 2021

# ClearSeal Screw-Cap Microtubes

- Ultra-clear polypropylene
- Packaged in resealable bags of 50
- Leak-proof O-ring seal
- Sterile, RNase, DNase & Pyrogen free
- White, printed graduations



C3205-C 0.5mL, sterile, w/O-ring & caps, non-graduated, conical bottom, 20 bags of 50 tubes, 1000/cs C3205-S 0.5mL, sterile, w/O-ring & caps, non-graduated, self-standing, 20 bags of 50 tubes, 1000/cs C3215-CG 1.5mL, sterile, w/O-ring & caps, non-graduated, conical bottom, 20 bags of 50 tubes, 1000/cs C3215-SG 1.5mL, sterile, w/O-ring & caps, printed graduations, self-standing, 20 bags of 50 tubes, 1000/cs 2.0ml, sterile, w/O-ring & caps, printed graduations, conical bottom, 20 bags of 50 tubes, 1000/cs C3220-CG C3220-SG 2.0mL, sterile, w/O-ring & caps, printed graduations, self-standing, 20 bags of 50 tubes, 1000/cs C3220-SGL \$239.50 2.0mL, sterile, w/O-ring & caps, printed graduations, self-standing, 20 bags of 50 tubes, 1000/cs

\$229.50 \$229.50 \$229.50 \$229.50 \$229.50 \$229.50



## **NEW!** Drawstring **Pull Closure**





two cases of bags (A8001R)



# BOWTIE Biohazard Bags

- Drawstring for easy closure
- 12-16 gallon capacity, 25 x 35in  $(45-60 \text{ liter capacity}, 64 \times 89 \text{ cm})$
- Color changing sterilization indicator
- Bin includes touchless lid with foot pedal





Sterilization Indicator

A8001R

BowTie™ Biohazard bags, HDPE, w/ drawstring, marking area and sterilization indicator, 100 bags/cs \$149.00

A8000B

BowTie™ Biohazard Bin, attached lid, w/ hands-free foot pedal, PP, 10 x 14 x 19in., 1/ea \*With purchase of two cases of BowTie™ bags (A8001R)

\$59.90 FREE!

3 inch

# FlipTop™ Freezer Boxes Cardboard and Polycarbonate

2 inch













All Boxes Include

• Standard  $5^{1}/4 \times 5^{1}/4 \times 2$  inch size

- Available in 6 colors and 81 or 100 place
- Polycarbonate models withstands 196°C
- Perfect for cryovials

#### 2 inch Boxes

Z IIICII DOXES					
2700	Cardboard freezer box, 2in., with dividers, 100-place, white, 5/pk				
2700-A*	Cardboard freezer box, 2in., with dividers, 100-place, rainbow, 5/pk				
2781	Cardboard freezer box, 2in., with dividers, 81-place, white, 5/pk				
2781-A*	Cardboard freezer box, 2in., with dividers, 81-place, rainbow, 5/pk				
2081	Polycarboante freezer box, 2in., 81-place, blue, 5/pk				
2100	Polycarbonate freezer box, 2in., 100-place, blue, 5/pk				
*For individual colors, replace (-A) with -B for Blue, -G for Green, -R for Red, -O for Orange, and -P for Purple					

Hinged lid

\$14.75

\$15.35 \$14.75 \$15.35 \$49.90 \$49.90



•  $5^{1}/4 \times 5^{1}/4 \times 3$  inch size

- Up to 12mm diameter tubes in 100-place
- Up to 13mm diameter tubes in 81-place
- 25-Place (5x5) for snap-cap 5mL MacroTubes<sup>®</sup>
- 36-Place (6x6) for screw-cap 5mL MacroTubes<sup>®</sup>

#### 3 inch Boxes

R3700	Cardboard freezer box, 3in., 100-place, white, 50/cs	\$139.00
R3781	Cardboard freezer box, 3in., 81-place, white, 50/cs	\$139.00
C2580-25	Cardboard freezer box, 3in., 25-place, for snap-cap 5mL tubes, 5/pk	\$19.80
C2580-36	Cardboard freezer box, 3in., 36-place, for screw-cap $5\mathrm{mL}$ tubes, $5/\mathrm{pk}$	\$19.80

## Alpha-numeric coding eliminates mix-ups


## nUVaClean UV pipette Carousel


- Simple, one-step decontamination
- Destroys >99% of stray DNA
- Store and decontaminate in one step
- Eliminates microbiological contaminants











BioWand™ Personal UV Sanitizer, 265-280nm, 1/ea \*With purchase of a nUVaClean (P5590)



# Halo<sup>™</sup> Pipettor Single Channel Adjustable Pipettes





Pipette

- 10-40% lighter than leading competitors
- 4 digit micrometer with DriftGuard™ protection
- Exclusive 1250µL for extended length tips
- User-enabled calibration port
- Compatible with universal tips

Cat. No.	Companible with universal tips	
H6800-2	Halo™ Single Channel Pipettor, lightweight, 0.2-2μL	\$249.00
H6800-10	Halo™ Single Channel Pipettor, lightweight, 1-10μL	\$249.00
H6800-20	Halo™ Single Channel Pipettor, lightweight, 2-20µL	\$249.00
H6800-50	Halo™ Single Channel Pipettor, lightweight, 5-50µL	\$249.00
H6800-100	Halo™ Single Channel Pipettor, lightweight, 10-100µL	\$249.00
H6800-200	Halo™ Single Channel Pipettor, lightweight, 20-200µL	\$249.00
H6800-1M	Halo™ Single Channel Pipettor, lightweight, 100-1000µL	\$249.00
H6800-1MXL	Halo™ Single Channel Pipettor, lightweight, 200-1250µL	\$249.00
H6800-5M	Halo™ Single Channel Pipettor, lightweight, 0.5-5mL	\$249.00
H6800-10M	Halo™ Single Channel Pipettor, lightweight, 1-10mL	\$249.00
H6800-CR	Carousel Stand for 6 Halo pipettes *With purchase of any 3 Halo Pipettors	\$109.50
	*With purchase of any 3 Halo Pipettors	FREE!*

# Microscope Slides Positively charged

• Positive charge - increased adhesion



Dispenser





M7100-45 Positive charge, 45° corners, frosted end, 25x75mm, 72/pk M7100-90 Positive charge, 90° corners, frosted end, 25x75mm, 72/pk

Microscope Slide Dispenser, holds 72 standard slides

\*With purchase of five packs of microscope slides (45° & 90°)

ProPette LE Pipettors
Multi-Channel Pipettors







- Wide volume range (1.0μL-300μL)
- Outstanding user comfort
- Compatible with standard tips
- Superior accuracy & reproducibility

P5208-10U 8-Channel Pipette, 1.0 to 10µL P5208-20U 8-Channel Pipette, 2.0 to 20µL \$399 P5208-50U 8-Channel Pipette, 5.0 to 50µL \$399 P5208-200U 8-Channel Pipette, 20 to 200µL \$399 P5208-300U 8-Channel Pipette, 30 to 300µL \$399 P5212-10U 12-Channel Pipette, 1.0 to 10µL \$499 P5212-20U 12-Channel Pipette, 2.0 to 20µL \$499 P5212-50U 12-Channel Pipette, 5.0 to 50µL \$499 P5212-200U 12-Channel Pipette, 20 to 200µL \$499 P5212-300U 12-Channel Pipette, 30 to 300µL P4406-P SureStand™ holds up to 4 multi-channels \$79 \*With purchase of any multi-channel pipette listed above

\$19.50

\$19.50

\$29.50

FREE!\*

Enhanced cell adhesion

### Pipette from 1μL to 300μL

# Centrifuge Tubes: 5 to 250mL

















C2530 5mL, non-sterile, screw-caps, 100 per bag, 500/cs C2535 5mL, sterile, screw-caps, 100 per bag, 500/cs C2540 5mL, sterile, screw-caps, foam racks of 50, 500/cs \$105.00 C2410 10mL, sterile, screw-caps, 20 bags of 50, 1000/cs \$195.00 C2600 15mL, sterile, screw-caps, foam racks of 50, 500/cs \$150.00 C2601 15mL, sterile, screw-caps, 10 bags of 50, 500/cs C2592 Workstation rack for 50 x 5, 10, & 15mL tubes, 1/ea \$39.50

With purchase of three cases of centrifuge tubes

- Vivid graduations
- A perfect size for any application
- Marking area on most sizes
- RNase, DNase, & Pyrogen free

• 25mL tube fits 50mL rotor Specifications:

All Tubes: -86°C to 121°C Temp. Range: Tubes: Polypropylene Caps: Polyethylene \*refers to liquids w/ a specific gravity of 1.0, using appropriate rotors, cushions and adapters.





C2625-B 25mL, sterile, screw-caps, 8 bags of 25, 200/cs C2625-R 25mL, sterile, screw-caps, foam racks of 25, 200/cs C2630 30mL, sterile, self-standing, 20 bags of 50, 500/cs C2602 50mL, sterile, screw-caps, foam racks of 50, 500/cs \$175.00 C2603 50mL, sterile, screw-cap, 20 bags of 50, 500/cs C2605 50mL, sterile, self-standing, 20 bags of 50, 500/cs \$179.00 C1250 250mL, sterile, screw-cap, 14 bags of 5, 70/cs \$195.00 R1070 Combo-rack for 30 x 15mL and 20 x 25/50mL, 5/pk \$59.00 \*With purchase of three case of centrifuge tubes

# NEW!

## DuoClick™ Screw-cap Culture Tubes with 2 position cap

Watch the Video!




	<u>12 x 75mm</u>	
T8730	5mL, polypropylene, sterile, 20 bags of 25 tubes, 500/cs	\$135.00
T8732	5mL, polypropylene, sterile, 10 foam racks of 50 tubes, 500/cs	\$159.00
T8734	5mL, polypropylene, sterile, individually wrapped, 500/cs	\$194.00
T8750	5mL, polypropylene, non-sterile,non-graduated, tubes and caps packed separately, 1000/cs	\$109.00


#### 17 x 100mm

T8830	14mL, polypropylene, sterile, 20 bags of 25 tubes, 500/cs	\$194.00
T8832	14mL, polypropylene, sterile, 10 foam racks of 50 tubes, 500/cs	\$239.00
T8834	14mL, polypropylene, sterile, individually wrapped, 500/cs	\$259.00
T8850	14mL, polypropylene, non-sterile, non-graduated, tubes and	\$209.00


- NEW DuoClick screw cap! (Patent Pending)
  - Aerobic culturing at first click
  - Anaerobic after second click
  - Easy to open Easy to close
- Standard Dual position caps available
- Sterile; RNase & DNase free







Cap is secure.



Cap is airtight

## 25mL Reservoirs & Lids | Microtube Racks

- Conserves even more reagent
- No-stick stacking design
- MicroTrough™ bottom
- 5 lids included in each case



Reservoir, 25mL, non-sterile, bulk packagaing, 100/cs P8025 \$52.00 P8025-1S Reservoir, 25mL, sterile, individually wrapped, 100/cs \$109.00 \$109.00 P8025-5S Reservoir, 25mL, sterile, 5 per sterile bag, 200/cs

NEW 25mL reservoir lids

P8025-LID 25mL solution reservoirs lids, sterile, 5 per sterile bag, 200/cs \$79.00

• 80 well (5 rows of 16) for 1.5 to 2.0mL microtubes

Rotate 1/4 turn

- Raised numbers on each well
- Recessed well bottoms for improved tube stability



Tube rack, fraction collector type, 80 x 1.5mL/2mL, rainbow 5/pk \$29.50 For individual colors, add -B for Blue, -G for Green, -R for Red, -O for Orange, and -Y for Yellow, to catalog number

# Transport Tubes



Shown in C2592 rack

- 95kPa rated
- Vivid, printed graduations
- Leak-proof, 1<sup>1</sup>/<sub>4</sub> turn screw-cap
   Sterile, with caps attached

C1811 5mL, screw-cap, 16 x 60mm, 100 per sterile bag, 1000/cs \$145.00 C1812 10mL, screw-cap, 16 x 80mm, 100 per sterile bag, 1000/cs \$169.00




C1812

CoolCaddy TM PCR Work Station



Watch the Video!





- Keeps samples cool, but not frozen
- Accomodates several different vessel types
  - Vials, Tubes, Strips, Plates & Reservoirs

Capacity: 1.5/2.0mL Microtubes: 2.0mL Cryotubes: 0.5mL Microtubes:

9 Tubes 8 Tubes

1 Plate 96 Tubes or 12 Strips

R4015

CoolCaddy™ PCR work station (17x 21 x 4.5cm)

C2592 Workstation rack for 50 tubes (16-17mm in diameter), 1/ea







