


# Round Bottom Snap Cap Tubes With Attached Caps

- Microtube convenience in larger capacities
- Attached caps seal at the rim with a “click”
- Sterile, easy-open bags of 100 tubes
- Molded graduations and frosted marking area



<b>MATERIAL:</b>	Polypropylene
<b>TEMPERATURE RANGE:</b>	-86°C - 121°C
<b>CENTRIFUGATION RATING:</b>	Up to 4,000×g
<b>DIMENSIONS:</b>	5mL - 12×55mm 7mL - 13×62mm 10mL - 14×82mm

## Stop fumbling with caps for your tubes!

MTC Bio’s new snap cap culture/centrifuge tubes adopt the conveniences of popular 1.5mL microcentrifuge tubes in higher capacities of 5mL, 7mL, and 10mL. Forget the cross-contamination concerns of mismatched caps!

These durable polypropylene tubes feature molded graduations on all sizes. The smaller 5mL and 7mL tubes also have a frosted marking area. Tubes are supplied in sterile, easy-open bags of 100 and are autoclavable.



CAT. NO.	DESCRIPTION	QTY	PRICE
<b>C0905</b>	5mL Sterile, w/ Attached Snap Cap, Molded Graduations, Frosted Marking Area	<b>1000</b>	<b>\$59.50</b>
<b>C0907</b>	7mL Sterile, w/ Attached Snap Cap, Molded Graduations, Frosted Marking Area	<b>1000</b>	<b>\$79.50</b>
<b>C0910</b>	10mL Sterile, w/ Attached Snap Cap, Molded Graduations	<b>1000</b>	<b>\$94.50</b>