

ProSeries™ Tubes

Performance Centrifuge Tubes

with
FreezeTech™
Temperature
Indicator



- 20,000×g max RCF marked directly on tube
- FreezeTech™ temperature indicator changes color at 0°C
- Hi-Viz™ graduations, including maximum freeze line
- Eco-Plastic Racks and/or resealable zip bags
- E-Beam sterilized (SAL 10⁻⁶)



20,000×g Maximum Centrifugation Rating



FreezeTech™ Indicator Changes Color at 0°C



Hi-Viz™ Graduations w/ Maximum Freeze Line



Available in Eco-Plastic Recyclable Racks

MTC Bio's new ProSeries™ Centrifuge Tubes offer maximum performance and a wide range of innovative features. ProSeries™ Tubes feature a 20,000×g RCF rating, an industry-first temperature indicator patch, and Hi-Viz™ Graduations that are recognizable at any viewing angle. Marked to both their maximum nominal capacity, as well as their maximum freeze volume, these tubes have been designed with a focus on environmentalism and waste reduction. Tubes are marked with standardized recycling information and are available in Eco-Plastic fold-flat racks rather than traditional styrofoam racks.

The flat, markable screw caps feature high precision threads that make them easy to open with one hand without compromising the 95kPa rated leak-proof seal.

CAT. NO.	DESCRIPTION	QTY	LIST PRICE
C2800	15mL, Sterile, with Temperature Indicator, 20 Eco-Plastic Racks of 25 Tubes	500	\$189.00
C2801	15mL, Sterile, with Temperature Indicator, 20 Bags of 25 Tubes	500	\$178.00
C2802	50mL, Sterile, with Temperature Indicator, 20 Eco-Plastic Racks of 25 Tubes	500	\$219.00
C2803	50mL, Sterile, with Temperature Indicator, 20 Bags of 25 Tubes	500	\$199.00