

ProSeries™ Tubes


High Performance Centrifuge Tubes

- 20,000xg Max RCF marked directly on tube
- FreezeTech™ indicator changes color at 0°C
- Hi-Viz™ graduations with max freeze line
- Resealable sterile bags and eco-plastic racks

Eco-Plastic Racks





Autoclavable, recyclable, and easily flattened after use
Also available for Standard 15 and 50mL Tubes - Page 2



15mL

50mL




A8403B **FREE BenchBin 3L™**
with any order of ProSeries™ Centrifuge Tubes - No Limit!
BenchBin™3L info on Page 4

Cat No.	Description	Qty	Price
C2800	15mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500	\$189.50
C2801	15mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500	\$178.00
C2802	50mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500	\$219.00
C2803	50mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500	\$199.00

Material
Tube: Polystyrene
Cap Frame: HDPE
Cap Mesh: Nylon

Flow Cytometry Tubes FACS Tubes and Cell Strainer Caps

- Compatible with FACS™ (BD) instruments
- Strainer Caps: 20µm, 35µm, 70µm, 100µm mesh
Standard, non-strainer caps also available
- 12x75mm (5mL) tubes in crystal clear polystyrene



Cat No.	Description	Qty	Price
Strainer and Standard Caps without Tubes			
T9009-20	20µm, Red, Strainer Cap, Translucent, Sterile, Bags of 25	500	\$200.25
T9009	35µm, Green, Strainer Cap, Translucent, Sterile, Bags of 25	500	\$200.25
T9009-70	70µm, White, Strainer Cap, Opaque, Sterile, Bags of 25	500	\$200.25
T9009-100	100µm, Yellow, Strainer Cap, Translucent, Sterile, Bags of 25	500	\$200.25
T9019	Standard Cap, Non-Strainer, Non-Sterile, Bags of 125	1000	\$82.80
FlowTubes™ - FACS Tubes with Included Caps			
T9005	5mL, 12x75mm, 35µm Strainer Cap, Sterile, Bags of 25	500	\$502.88
T9005-1S	5mL, 12x75mm, 35µm Strainer Cap, Sterile, Individual Wrap	500	\$571.32
T9010	5mL, 12x75mm, Standard Cap, Sterile, Bags of 25	500	\$118.64
FlowTubes™ - FACS Tubes without Caps			
T9014	5mL, 12x75mm, Without Cap, Sterile, Tray Packs of 250	1000	\$144.48
T9020	5mL, 12x75mm, Without Cap, Non-Sterile, Bulk Econo-Pack	4000	\$259.00

Flip 'N Fold™

Sterile Sample Bags

- Quick and easy to open and reseal your samples
- Leak-proof and useable between -210°C and 80°C
- Virgin polyethylene for maximum strength and clarity
- Sterilize RNase, DNase, and Pyrogen free

Marking areas on select sizes



55oz - 300×180mm



2oz - 127×76mm



4oz - 178×76mm



18oz - 229×114mm



27oz - 304×114mm



36oz - 381×114mm

Cat No.	Description	Qty	Price
B5750	2oz Sampling Bag, 127mm×76mm, Sterile, 2.8mil Thickness	500	\$80.40
B5760	4oz Sampling Bag, 178mm×76mm, Sterile, 2.8mil Thickness	500	\$83.40
B5761	4oz Sampling Bag, 127mm×76mm, Sterile, 2.8mil Thickness, Marking Area	500	\$90.00
B5770	18oz Sampling Bag, 229mm×114mm, Sterile, 2.8mil Thickness	500	\$109.80
B5771	18oz Sampling Bag, 229mm×114mm, Sterile, 2.8mil Thickness, Marking Area	500	\$138.00
B5780	27oz Sampling Bag, 304mm×114mm, Sterile, 2.8mil Thickness	500	\$166.80
B5790	36oz Sampling Bag, 381mm×114mm, Sterile, 2.8mil Thickness	500	\$186.00
B5910	55oz Sampling Bag, 300mm×180mm, Sterile, 4.0mil Thickness	250	\$178.80
B5911	55oz Sampling Bag, 300mm×180mm, Sterile, 4.0mil Thickness, Marking Area	250	\$210.00

ClearSeal™

Screw Cap Microtubes

- Leak-proof O-Ring seal
- Ultra-clear polypropylene
- Resealable sterile bags of 50
- RNase, DNase, and Pyrogen free
- White printed graduations
- Conical or standing bottoms



0.5 mL Conical



0.5 mL Self-Standing



1.5 mL Conical



1.5 mL Self-Standing



2.0 mL Conical



2.0 mL Self-Standing



2.0 mL Self-Standing Loop Cap

Cat No.	Description	Qty	Price
C3205-C	0.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000	\$247.57
C3205-S	0.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap,	1000	\$247.57
C3215-CG	1.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000	\$247.57
C3215-SG	1.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000	\$247.57
C3220-CG	2.0mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap, Graduated	1000	\$247.57
C3220-SG	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000	\$247.57
C3220-SGL	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Loop Screw Cap, Graduated	1000	\$256.27

Temperature
Min: -86°C
Max: 121°C

Material
Polypropylene

Max RCF
22,000×g

Centrifuge Tubes

5mL to 500mL Capacities


- Vivid graduations
- Perfect size for any application
- RNase, DNase, & pyrogen free

Temperature
Min: -86°C
Max: 121°C

Material
Tube: Polypropylene
Cap: Polyethylene



FREE R1030 Rack
with purchase of any tube in this group



Cat No.	Description	Qty	Price
5mL MacroTubes®			
C2530	Separate Caps, Non-Sterile, Bulk Bagged	500	\$99.00
C2535	Attached Caps, Sterile, Bags of 100	500	\$105.00
C2540	Attached Caps, Sterile, Racks of 50	500	\$112.35
C2540-OB	Attached Caps, Opaque Black, Sterile, Racks of 50	500	\$129.50
C2545	Attached Caps, Vivid Graduations, Sterile, Racks of 50	500	\$119.00
10mL Centrifuge Tubes			
C2410	Attached Cap, Sterile, Bags of 50	1000	\$234.00
15mL Centrifuge Tubes			
C2600	Attached Cap, Sterile, Racks of 50	500	\$172.50
C2601	Attached Cap, Sterile, Bags of 50	500	\$162.00
C2607	Attached Cap, Sterile, Eco-Plastic Racks of 25	500	\$172.50
C2615-W	Attached Cap, Sterile, Individually wrapped	500	\$117.47

Cat No.	Description	Qty	Price
25mL Centrifuge Tubes			
C2625-B	Attached Cap, Sterile, Bags of 25	200	\$101.20
C2625-R	Attached Cap, Sterile, Racks of 25	200	\$104.65
30mL Centrifuge Tubes			
C2630	Screw Cap, Self-Standing, Sterile, Bags of 25	500	\$169.05
50mL Centrifuge Tubes			
C2602	Attached Cap, Sterile, Racks of 25	500	\$199.50
C2603	Attached Cap, Sterile, Bags of 25	500	\$179.40
C2603-B	Separate Screw Cap, Non-Sterile, Bulk Bag	500	\$164.34
C2608	Attached Cap, Sterile, Eco-Plastic Racks of 25	500	\$199.50
C2650-W	Attached Cap, Sterile, Individually Wrapped	300	\$143.06
C2605	Attached Cap, Sterile, Self-Standing, Bags of 25	500	\$205.85

Cat No.	Description	Qty	Price
Large Volume Centrifuge Tubes			
C0175	175mL, Sterile PP, Conical Bottom, Bags of 6	120	\$232.00
C0225	225mL, Sterile, PP, Conical Bottom, Bags of 5	90	\$195.00
C1250	250mL, Sterile, PP, Conical Bottom, Bags of 5	70	\$234.00
C1500	500mL, Sterile, PP, Conical Bottom, Bags of 5	40	\$239.00
Centrifuge Tube Racks			
R1030	4×50/25mL, 12×15/10/5mL, 32×1.5/0.5mL Tubes	5	\$63.13
C2592	Decanting, 50×5mL, Silicone Grip Wells	4	\$47.40
R1070	30 Slots for 10/15mL, 20 Slots for 25/50mL	5	\$63.13
R1900	6 Slots for 175, 225, or 250mL Tubes	5	\$69.00
C1500-RK	4 Slots for 500mL Tubes	1	\$86.11

Solution Reservoirs

Standard 12-Channel Width

- Graduated inner walls
- Micro-Trough U-shaped bottoms
- Pouring spouts on all 4 corners
- RNase, DNase, & pyrogen free
- Pack of 5 protective lids with standard 25 and 50mL reservoirs



Description	25mL* Standard	50mL* Standard	50mL Divided	100mL Standard	QTY	Price
Non-Sterile Bulk Pack	P8025	P8050	P8051	P8100	100	\$58.24
Sterile Individually Wrapped	P8025-1S	P8050-1S	P8051-1S	P8100-1S	100	\$122.08
Sterile Packs of 5	P8025-5S	P8050-5S	P8051-5S	P8100-5S	200	\$122.08
Lids - Sterile Packs of 5	P8025-LID	P8050-LID	N/A	N/A	200	\$84.53

*5 FREE Lids included with every purchase of Standard 25mL and Standard 50mL reservoirs.

ASPIR-8™ Reservoirs

Reduced 8-Channel Width

- Specifically designed for 8-channel pipettes
- 75mm inner chamber
- Conserve valuable reagents



Description	5mL 8-Channel	10mL 8-Channel	25mL 8-Channel	QTY	Price
Non-Sterile Bulk Pack	P7005	P8010	P7025*	300*	\$111.28*
Sterile Individually Wrapped	P7005-1S	P8010-1S	P7025-1S	100	\$111.28
Sterile Packs of 5	P7005-5S	P8010-5S	P7025-5S	200	\$111.28

*P7025 has quantity of 100/pack and price of \$55.64

Sep-R-8™ Plates

Break-Away PCR Plates

- Breaks into 12 strips of 8 PCR tubes
- Ultra-thin walled tube wells
- Alpha-numeric reference grids

Capacity
96 × 0.2mL

Material
Polypropylene

Style
Non-Skirted




FREE Tool (#P0800-TL) 2-Pack with orders of 2 packs of Sep-R-8™

Cat No.	Description	Qty	Price
P9602-N-T	Sep-R-8™ Break-Away PCR Plate, 96×0.2mL, Non-Skirted, Breaks into Strips of 8 PCR Tubes	50	\$149.00
P0800-TL	Tool for Opening/Closing PCR Tubes and Strips, with Flat Edge for Sealing Plates	2	\$86.80

OptiWell™ Deep Well Plates


24, 48, 96, and 384 Well Plate Options

- For automated research platforms and robotic equipment
- Medical-grade polypropylene - minimal binding / liquid retention
- Black printed alpha-nums on 24, 48, and 96 well plates
- Choose from round (U), conical (V), or flat (F) bottoms



1 + 1 = 3 Buy 2 Get 1 FREE
Same item only

Scan this QR code for over 50 Deep Well Plate configurations, sealing mats, and combs



Cat No.	Description	Qty	Price
D3024-01S	24 Square Well, 10mL, 41.9mm Height, V-Bottom, Sterile Bags of 5	50	\$258.00
D3048-01S	48 Square Well, 4.6mL, 41.2mm Height, U-Bottom, Sterile Bags of 5	50	\$280.00
D3096-01S	96 Round Well, 0.36mL, 14.75mm Height, V-Bottom, Sterile Bags of 10	100	\$520.00
D3096-01S	96 Round Well, 0.4mL, 14.75mm Height, F-Bottom, Sterile Bags of 10	100	\$520.00
D3096-09S	96 Round Well, 0.5mL, 15.1mm Height, V-Bottom, Sterile Bags of 5	50	\$200.00
D3096-11S	96 Round Well, 0.5mL, 14.5mm Height, U-Bottom, Sterile Bags of 10	100	\$520.00
D3096-05S	96 Round Well, 1.0mL, 33.0mm Height, U-Bottom, Sterile Bags of 5	50	\$276.00
D3096-25S	96 Square Well, 1.2mL, 24.0mm Height, U-Bottom, Sterile Bags of 5	50	\$200.00
D3096-39S	96 Square Well, 2.2mL, 44.2mm Height, V-Bottom, Sterile Bags of 5	50	\$200.00
D3096-41S	96 Square Well, 2.2mL, 44.2mm Height, U-Bottom, Sterile Bags of 5	50	\$200.00
D3384-01S	384 Square Well, 120µL, 14.6mm Height, V-Bottom, Sterile Bags of 10	100	\$676.00

Promo Code: RE2024

Hinged-Lid Freezer Boxes

Temperature
Min: -196°C

Footprint
5¼ × 5¼ inches

Polycarbonate or Cardboard Boxes

- Attached hinged lids
- Alpha-numeric reference grids
- 9×9 or 10×10 formats
- LN₂ drainage slots



3" Cardboard
81-Place



2" Polycarbonate
100-Place



2" Cardboard
100-Place


Cat. No.	Description	Qty	Price
Cardboard Freezer Boxes - 5¼×5¼ Footprint, Attached Hinged Lid			
R2700*	2" Height, 100-Place, Alphanumeric Grids, White Color*	5	\$18.59
R2781*	2" Height, 81-Place, Alphanumeric Grids, White Color*	5	\$18.59
R3700	3" Height, 100-Place, Alphanumeric Grids, White Color	50	\$148.73
R3781	3" Height, 81-Place, Alphanumeric Grids, White Color	50	\$148.73
Polycarbonate Freezer Boxes - 5¼×5¼ Footprint, Attached Hinged Lid			
R2200	2" Height, 100-Place, Alphanumeric Grids, Transparent Lid	5	\$59.88
R2281	2" Height, 81-Place, Alphanumeric Grids, Transparent Lid	5	\$59.88

*Add "-B" for Blue, "-G" for Green, "-O" for Orange, "-P" for Purple, "-R" for Red, or "-AST" for Assorted Colors. Price: \$19.34

Round Bottom Snap Cap Tubes

Hinged Cap Centrifuge Tubes

- Microtube convenience in larger capacities
- Attached caps seal at the rim with a "click"
- Sterile, easy-open bags of 100 tubes



Temperature
Min: -86°C
Max: 121°C

Material
Polypropylene

Max RCF
4000×g

Cat. No.	Description	Qty	Price
C0905	5mL, Attached Snap Cap, Molded Graduations, Sterile Bags of 100	1000	\$59.50
C0907	7mL, Attached Snap Cap, Molded Graduations, Sterile Bags of 100	1000	\$79.50
C0910	10mL, Attached Snap Cap, Molded Graduations, Sterile Bags of 100	1000	\$94.50

BenchBin™ Biohazard Bin

Benchtop Wastebin with Spring Lid

- Secures and covers hazardous waste
- Safely dispose of tips, tubes, gloves

FREE
BenchBin 3L™
Details on page 1

A9004R
8½×11" Bags

A9006R
10×12½" Bags



1.5L Capacity



3L Capacity

Cat. No.	Description	Qty	Price
A8401B	BenchBin™, 7×5×7", 1.5L Capacity, fits A9004R Autoclave Bags	1	\$39.80
A8403B	BenchBin 3L™, 6×5×8½", 3L Capacity, fits A9006R Autoclave Bags	1	\$47.50
A9004R	Autoclave Bags, 8½×11" Printed Markings, fits A8401B BenchBin™	100	\$31.00
A9006R	Autoclave Bags, 10×12½" Printed Markings, fits A8403B BenchBin 3L™	100	\$39.00
A8401-SK	Starter Kit: 1 BenchBin™ 1.5L (A8401B) and 400 Autoclave Bags (A9004R)	1	\$120.75

Erlenmeyer Flasks

Flat or Baffled Bottoms

- Polycarbonate or PETG
- Vented Screw caps w/ 0.2µL filter
- Sterile, individually wrapped
- RNase, DNase, and pyrogen free




Buy 2 Get 1 FREE Free flask must be of equal or lesser value

Cat. No.	Description	Qty	Price
Erlenmeyer Shaker Flasks w/ Vented Screw Caps - PETG, Flat Bottom			
F4060-F	125mL, PETG, Flat Bottom, Sterile, Individually Wrapped,	24	\$213.60
F4061-F	250mL, PETG, Flat Bottom, Sterile, Individually Wrapped,	12	\$120.00
F4062-F	500mL, PETG, Flat Bottom, Sterile, Individually Wrapped,	12	\$154.80
F4063-F	1000mL, PETG, Flat Bottom, Sterile, Individually Wrapped,	6	\$117.00
Erlenmeyer Shaker Flasks w/ Vented Screw Caps - PC, Baffled Bottom			
F4060-B	125mL, PC, Baffled Bottom, Sterile, Individually Wrapped,	24	\$239.00
F4061-B	250mL, PC, Baffled Bottom, Sterile, Individually Wrapped,	12	\$136.00
F4062-B	500mL, PC, Baffled Bottom, Sterile, Individually Wrapped,	12	\$179.00
F4063-B	1000mL, PC, Baffled Bottom, Sterile, Individually Wrapped,	6	\$124.40
Jumbo Fernbach Flasks w/ Vented Screw Caps - Styrene-Butadiene, Baffled Bottom			
F4086	2.8L, Styrene-Butadiene, Baffled Bottom, Sterile, Individually Wrapped	6	\$270.00
F4087	5.0L, Styrene-Butadiene, Baffled Bottom, Sterile, Individually Wrapped	4	\$380.00


Microscope Slides

Positive Charged for Cell Adhesion




90° Corners

- Choose from 90° or 45° corners
- Ultra-clean / dust-free
- Frosted on one end
- Dispenser fits 72 slides



FREE
Dispenser
M7100-SD
with orders of
5 slide packs



45° Corners

Cat. No.	Description	Qty	Price
M7100-90	90° Corners, Positive Charged, Frosted End, 25×75mm	72	\$21.84
M7100-45	45° Corners, Positive Charged, Frosted End, 25×75mm	72	\$21.84
M7100-SD	Dispenser for 72 Standard 25×75mm Microscope Slides	1	\$38.23

BowTie™ Biohazard Family

Drawstring Bags and Hands-Free Bins

- Easy drawstring closure
No wire tie needed!
- 12 or 30 gallon capacity bags
- Color-changing sterilization indicator
- Hands-free foot pedal lid on bins



White Marking Area
Sterilization Indicator
Before Autoclave
After Autoclave
Biohazard Symbol



30 Gallon
33×41" Bag



26 Gallon
19×19×33" Bin

Cat. No.	Description	Qty	Price
BowTie™ Biohazard Bags - Drawstring, Marking Area, Sterilization Indicator			
A8001R	12 Gallon, 25×35", Fits A8000B Biohazard Bin	100	\$159.43
A8003R	30 Gallon, 33×41", Fits A8002B Biohazard Bin	100	\$187.50
BowTie™ Biohazard Bins - Attached Lid and Hands-Free Foot Pedal			
A8000B	8 Gallon, 10×14×19", Fits A8001R Biohazard Bag	1	\$78.53
A8002B	26 Gallon, 19×19×33", Fits A8003R Biohazard Bag	1	\$189.00



Bags are oversized for a perfect fit when folded over the rim

To receive promotional packages, please reference Promo Code RE2024 on your order. All promotions are valid while supplies last and expire December 31, 2024. Prices are subject to change without notice.

Your Personalized Contact Footer Here!